

Assay of Complex IV (cytochrome C oxidase)

Principle:

- Complex IV in the respiratory chain transfers electrons from reduced cytochrome c to oxygen. Its activity will be assessed by following the decrease of reduced cytochrome c absorbance at 550 nm.

Practical set up:

- 6 samples may be analyzed on the same day (5 patients and 1 control).
- Reaction medium composition:
 - 50 mM K phosphate pH 7.0 (25 mM for skin fibroblasts)
 - 100 μ M reduced cytochrome c (50 μ M for skin fibroblasts)
 - Tissue: 40 μ g of proteins (post-nuclear supernatant (liver or muscle) or cell suspension) or 1 μ g of mitochondrial proteins
- Preparation of reaction medium:
 - 1) In a 15 mL tube, prepare 12 to 13 mL of 100 μ M cytochrome c in 50 mM K phosphate pH 7.0
 - 2) Prepare a “**100% oxidized solution**” with 1mL of initial cytochrome c solution oxidized with few grains of **potassium ferricyanide** in a 1 mL cuvette (the color of the solution turns to a brownish dark red).
 - 3) Prepare a “**100% reduced solution**” with 1mL of initial cytochrome c solution reduced with few grains of **sodium dithionite** in a 1 mL cuvette (the color of the solution turns to salmon pink).
 - 4) Read the absorbance of the “100% oxidized solution” at 550 nm in the spectrophotometer after having obtained a blank on air. The absorbance should be around 0.7 for the 100 μ M solution.
 - 5) Re-blank the spectrophotometer on the “100% oxidized solution” and read the absorbance of the “100% reduced solution”. Its value is considered 100% reduction.
 - 6) Transfer an aliquot (50 or 100 μ L) of the “100% reduced solution” in the initial cytochrome c solution and read its new absorbance. Keep on adding aliquots until reaching an absorbance between 90 and 95% of the 100% reduced solution absorbance. For example, if the 100% reduced solution had an absorbance of 1.4, the initial cytochrome c solution must be progressively reduced until its absorbance is between 1.26 and 1.33.
- Assay:
 - 1) Reading in the spectrophotometer, at 37°C, at wavelength 550 nm
Initial calibration is performed on air.
 - 2) Incubate the cuvettes containing 980 μ L of reduced initial solution of cytochrome c **at 37°C, during 5 min**, in the spectrophotometer
 - 3) Start the reaction by adding **20 μ L of supernatant or cell suspension** (with proteins concentration set at 2 mg/mL) or of **beef heart mitochondria** (diluted at 0.05 mg/mL i.e. 1/800).
 - 4) Reading every 10 seconds during 3 minutes,
Two cuvettes may be read at the same time.
If the decrease appears too rapid giving a non-linear curve, re-do the assay with less tissue.
- Calculation:
 - 1) Complex IV specific activity is calculated as nanomoles/min/mg of post-nuclear supernatant or cell suspension proteins.
 - 2) The extinction coefficient for cytochrome c is $\epsilon = 18.5$
 - 3) The correction factor is therefore 1351.4 for 40 μ g of proteins in the assay (tissue homogenate) and 54054 for 1 μ g of proteins in the assay (isolated mitochondria)