ATELIER

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### **Recent publications:**

- Candat A, Paszkiewicz G, Neveu M, Gautier R, Logan DC, Avelange-Macherel MH, Macherel D (2014) The ubiquitous distribution of LEA proteins across cell compartments offers tailored protection against abiotic stress. *Plant Cell*, in press.
- Benamar A, Pierart A, Baecker V, Avelange-Macherel MH, Rolland A, Gaudichon S, di Gioia L, Macherel D (2013) A simple system using
  natural mineral water for high-throughput phenotyping of Arabidopsis thaliana seedlings in liquid culture. Int. J. H-T Screening. 4: 1-15.
- Candat A, Poupart P, Andrieu JP, Chevrollier A, Reynier P, Rogniaux H, Avelange-Macherel MH, **Macherel D** (2013) Experimental determination of organellar targeting-peptide cleavage sites using transient expression of green fluorescent protein translational fusions. *Anal. Biochem.* 434: 44-51.

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My first degree was from the University of Dundee, Scotland, followed by a PhD at University College London. I then had postdoctoral positions in Kent, England and Montpellier, France, in cell physiology and molecular biology before starting research into mitochondrial cell biology as a postdoctoral fellow at Oxford University in 1996. In Oxford I was frustrated with not being able to see mitochondria unambiguously in living tissue so generated stable transgenics of the model plant Arabidopsis expressing mitochondrial targeted GFP. In 1998 I left academic research and worked briefly for the UK government in London before returning to mitochondrial research in 2000 at the University of St Andrews, Scotland. In St Andrews I mutagenized my mito-GFP line and screened, by epifluoresence microscopy, for mutants with altered mitochondrial shape, size, number or cellular distribution. In 2008 I moved to Canada as Associate Professor before coming to Angers in 2012.

### **Recent publications:**

- Candat A, Paszkiewicz G, Neveu M, Gautier R, Logan DC, Avelange-Macherel MH, Macherel D (2014) The ubiquitous distribution of LEA
  proteins across cell compartments offers tailored protection against abiotic stress. *Plant Cell*, in press.
- Ren M, Venglat P, Qiu S, Feng Li, Cao Y, Wang E, Xiang D, Wang J, Alexander D, Chalivendra S, Logan DC, Matoo A, Gopalan S, Datla R (2012) Target of rapamycin signaling regulates metabolism, growth and lifespan in Arabidopsis. *Plant Cell* 24: 4850 4874.
- Schwarzländer M, Logan DC, Johnston IG, Jones NS, Meyer AJ, Fricker MD, Sweetlove LJ (2012) Pulsing of membrane potential in individual mitochondria: A stress-induced mechanism to regulate respiratory bioenergetics in Arabidopsis. Plant Cell 24: 1188-1201.

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### **Stephan WAGNER**

### PostDoc in Plant Energy Biology

Chemical Signalling Group, Institute for Crop Sciences and Ressource conservation (INRES), University of Bonn (Germany)

Homepage: http://www.plant-energy.uni-bonn.de/ Contact: <a href="mailto:stephan.wagner@uni-bonn.de">stephan.wagner@uni-bonn.de</a>

Using genetically encoded biosensors, we seek to understand how plants as sessile organisms react to a changing environment. Using a combination of *in vitro*, *in vivo* and studies on isolated mitochondria, we use our strong expertise in bioimaging to analyse how specifically redox and calcium signals shape mitochondrial respiration and metabolism. Before joining Markus Schwarzländer's lab in Bonn, I solved three-dimensional structures of proteins that are essential for a plant's resistance against pathogens. I am thus deeply interested in how these signals are mechanistically transduced to and from mitochondria and whether they regulate key mitochondrial proteins through direct modifications.

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### 5<sup>e</sup> ATELIER DU RESEAU MEETOCHONDRIE

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# Investigating plant mitochondrial dynamics by wide-field fluorescence microscopy

### David C Logan & David Macherel

The purpose of the plant side to the workshop is to compare and contrast confocal and widefield microscopy of plant mitochondria but using different experimental materials. Whilst the confocal part will focus on the use of mitochondrial targeted photo-convertible mEos (Mathur *et al*, 2010) to investigate inter-mitochondrial interactions, the wide-field component will focus on mitochondria and the actin cytoskeleton, and separately the visualisation of mitochondrial DNA (stained with SYBR Green I) and mitochondria (stable transgenics expressing mitochondrial targeted mCherry fluorescence protein.

Cytoskeleton and mitochondria visualization To visualize actin filaments and mitochondria, a double transgenic line was generated by crossing mito-GFP (Logan & Leaver, 2000) with a second line (mCherry-mTn) expressing an inframe fusion of mCherry (Shaner et al, 2004) to the C-terminal 197 amino acids of mouse talin which forms the actin binding domain (ABD). The mCherry-mTn line was generated by transformation of Arabidopsis Col-0 with pDCLmcherry-mTn. The backbone of pDCLmcherrymTn, the pDCLmcherry-X destination vector, was a modified pMDC43 (Curtis & Grossniklaus, 2003) in which mGFP6 was replaced by mCherry via pMDC7 and pMDC24. A cDNA fragment encoding the mTalin ABD was PCR-amplified from the pZP202-GFP-mTn vector (see Kost et al. 1998) using attB Gateway primers and recombined into pDONRz. Recombination between pDONRz-mTn and pDCLmCherry-X created pDCLmcherry-mTn. Actin and mitochondria will be imaged using a Zeiss Axioimager Z1.M with Colibri LEDs and a Hamamatsu Orca Flash 4.0 v2 camera. The system is controlled by Zeiss Zen software running on a Dell workstation.

# mtDNAandmitochondrialdynamicsA stable transgenic Arabidopsis line expressing mito-mCherry will be used (Candat *et al.* ThePlant Cell, in press). Mitochondrial DNA in this line will be stained using SYBR Green I(Invitrogen). Seedlings will be placed in 1 ml of MS culture medium and 0.2 μL of SYBR Greenstock will be added (i.e. 1:5000). Incubate for 15 min, wash 2x.

### Photoconvertible GFP as tool for mitochondrial dynamics

We will use a stable transgenic line of Arabidopsis constitutively expressing mitochondrialtargeted mEosFP (Mathur *et al.*, 2010). EosFP is a coral FP that was genetically engineered to yield a monomeric green-to-red photoconverting FP (Wiendemann *et al.*, 2004). Seedlings grown for 7-14 days on 0.5x MS agar medium will be observed using a Nikon confocal A1 system (GFP 488/525, RFP 561/595) and conversion will be induced with a 405 nm laser.

# Visualisation of mitochondrial membrane potential with TMRM in Arabidopsis (optional – time permitting)

Transgenic Arabidopsis seedlings constitutively expressing GFP in mitochondria (Logan *et al.* 2000; Candat *et al.*, 2014) will be grown in liquid culture for 10 days, and incubated with 20 nM TMRM [Tetramethylrhodamine methyl ester] for 15-30 min. Mitochondria labeled with GFP/TMRM will be observed with the Nikon confocal A1 system (GFP 488/525, RFP 561/595).

- Candat A, Paszkiewicz G, Neveu M, Gautier R, Logan DC, Avelange-Macherel MH, Macherel D (2014) The ubiquitous distribution of LEA proteins across cell compartments in Arabidopsis offers tailored protection against abiotic stress. Plant Cell, in press.
- Curtis MD & Grossniklaus U (2003) A gateway cloning vector set for high-throughput functional analysis of genes in planta. *Plant Physiol* **133**: 462–469
- Kost B, Spielhofer P & Chua NH (1998) A GFP-mouse talin fusion protein labels plant actin filaments in vivo and visualizes the actin cytoskeleton in growing pollen tubes. *Plant J* **16:** 393–401
- Logan DC & Leaver CJ (2000) Mitochondria-targeted GFP highlights the heterogeneity of mitochondrial shape, size and movement within living plant cells. *J Exp Bot* **51**: 865–871
- Mathur J, Radhamony R, Sinclair AM, Donoso A, Dunn N, Roach E, Radford D, Mohaghegh PSM, Logan DC, Kokolic K & Mathur N (2010) mEosFP-based green-to-red photoconvertible subcellular probes for plants. *Plant Physiol* **154**: 1573–1587
- Shaner NC, Campbell RE, Steinbach PA, Giepmans BNG, Palmer AE & Tsien RY (2004) Improved monomeric red, orange and yellow fluorescent proteins derived from Discosoma sp. red fluorescent protein. *Nat Biotechnol* 22: 1567–1572
- Wiendenmann J, Ivanchenko S, Oswald F, Schmitt F, Rocker C, Salih A, Sîndler KD, Nenhaus GU (2005) EosFP, a fluorescent marker protein with UV-inducible green-to-red fluorescnce conversion. Proc Natl Acad Sci USA **101**: 15905-15910

### 5° ATELIER DU RESEAU MEETOCHONDRIE

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### Image presentation and analysis with Fiji (ImageJ)

### David C Logan

It is recommended to use the latest version of Fiji. The acronym Fiji stands for "Fiji Is Just ImageJ" and is a distribution package of ImageJ complete with many plugins (http://fiji.sc/Fiji).

### Contents

- 1. Opening images.
  - 1.1. Bio-Formats Import Options.
- 2. Preparation of publication/presentation ready images.
  - 2.1. Adjustment of "Brightness and Contrast".
  - 2.2. Adding colour to a single grayscale image.
  - 2.3. Working with channels creating a merged image.
  - 2.4. Cropping an image.
  - 2.5. Addition of a scale bar.
  - 2.6. Addition of a time stamp.
  - 2.7. Addition of a z-position stamp.
  - 2.8. Adding annotations to the image using ROI manager.
  - 2.9. Exporting as a movie file
- 3. Analysis of mitochondrial dynamics qualitative and quantitative.
  - 3.1. Co-localisation for qualitative and quantitative representation of dynamics.
  - 3.2. Counting and measuring particles (mitochondria).

### 1. Opening images.

On a Mac, assuming the file suffix has been correctly attributed to Fiji then double clicking on the file icon with open the file.

On a PC, such as in the computer room at the workshop, select "File" then "Open" and then select the correct file. Alternatively, select the file icon, right-click, and choose "Open with" "ImageJ".

On up to date versions of Fiji different types of file will be opened automatically. However, with older versions of Fiji/ImageJ, Nikon .nd2 files (as saved automatically by the Elements software running the confocal) cannot be opened. Even with the newest version .nd2 files must be viewed as "Standard ImageJ" files (select this option in the "Bio-Formats Import Options" dialogue box, also, make sure to check "Split Channels".

### 1.1. Bio-Formats Import Options.

This dialogue box allows you to select various parameters to apply to the image when it is opened in Fiji. The newest version of Fiji (currently 1.49b) can handle most image files as a "Hyperstack" and this is the default setting. Typically, other selections to make are "Autoscale" (bottom left of dialogue) and "Split Channels", however, the choice of these and other selections will depend on the precise analysis to be performed.

Stack viewing		Metadata viewing	Information
View stack with:	Hyperstack ‡	🗆 Display metadata	Concatenate series when compatible - Allows multiple image series to be joined end to end.
Stack order:	XYCZT ‡	Display OME-XML metadata	Example: You want to join two sequential
		Display ROIs	timelapse series.
Dataset organiza	ation	Memory management	
Group files wit	h similar names	Use virtual stack	
Open files indi	vidually	Specify range for each series	
Swap dimensions		Crop on import	
Open all series	;		
Concatenate se	eries when compatible	Salit into separate windows	
Stitch tiles		Split channels	
Color options		Split focal planes	
Color mode:	Default \$	Split timepoints	
Autoscale			

"Autoscale" ensures that the histogram for each channel is stretched to each channels min and max values throughout the stack – adjustments can still be made using "Brightness and Contrast", but selecting "Autoscale" is a sensible start.

"Split Channels" is useful since often one wants to make a merged image and apply particular LUTs to the separate channels that are to be merged. A LUT (Look Up Table) is simply a colour-map array (palette) used to specify the colour and intensity of the image. After opening the image one is ready to proceed with preparation of a publication/presentation-ready copy, or some form of image analysis. Each of these will be discussed in turn.

### 2. Preparation of publication/presentation ready images.

<u>All modifications to an image should be performed on a copy and the original file</u> <u>direct from the image capture software should be saved permanently unaltered</u>. Many changes to images in Fiji cannot be undone and therefore at each stage when a new annotation or adjustment is to be applied it is recommended that a copy of the image or stack is made ("Image", "Duplicate"), make sure the "duplicate stack" click box is selected if duplicating a stack of images. Preparation of publication/presentation ready copies may involve adjustment to the "Brightness and Contrast" of the image, applying colour to a grayscale image, addition of a scale bar, time stamp, or z-position stamp, and annotation of objects within the image.

### 2.1. Adjustment of "Brightness and Contrast".

The following text is copied from the ImageJ guide:

"Use this tool to interactively alter the brightness and contrast of the active image. With 8bit images, brightness and contrast are changed by updating the image's lookup table (LUT), so pixel values are unchanged. With 16-bit and 32-bit images, the display is updated by changing the mapping from pixel values to 8-bit display values, so pixel pixel values are also unchanged. Brightness and contrast of RGB images are changed by modifying the pixel values.

**Histogram** The line graph at the top of the window, which is superimposed on the image's histogram, shows how pixel values are mapped to 8-bit (0--255) display values. The two numbers under the plot are the minimum and maximum displayed pixel values. These two values define the display range, or 'window'. ImageJ displays images by linearly mapping pixel values in the display range to display values in the range 0--255. Pixels with a value less than the minimum are displayed as black and those with a value greater than the maximum are displayed as white.

*Minimum and Maximum sliders* Control the lower and upper limits of the display range. Holding down Shift will simultaneously adjust all channels of a composite image (e.g., File⊳Open Samples⊳HeLa Cells (1.3 M, 48-bit RGB)).

**Brightness slider** Increases or decreases image brightness by moving the display range. Holding down Shift will simultaneously adjust all

channels of a composite image.

*Contrast slider* Increases or decreases contrast by varying the width of the display range.

The narrower the display range, the higher the contrast. Holding down Shift will simultaneously adjust all channels of a composite image.

**Auto** ImageJ will automatically optimize brightness and contrast based on an analysis of the image's histogram. Create a selection, and the entire image will be optimized based on an analysis of the selection. The optimization is done by allowing a small percentage of pixels in the image to become saturated (displayed as black or white). Each additional click on *Auto* increases the number of saturated pixels and thus the amount of optimization.

**Reset** Restores the original brightness and contrast settings. The display range is set to the full pixel value range of the image. A resetMinAndMax() macro call is generated if the command recorder is running. Holding down Shift restores original settings in all channels of a composite image.

*Set* Allows to enter the minimum and maximum display range values in a dialog box. A setMinAndMax() macro call is generated if the command recorder command recorder is running.

A 16-bit image consists of 65536 possible gray levels. Most of times, however, the relevant image information is contained only within a narrow range of the grayscale. This is the case, e.g., in low light microscopy, in which signal is restricted to the lower end of the grayscale. The *Set Display Range* dialog allows you to choose how to scale the range of gray levels of 16-bit images."

### N.B.

(a) A note of caution: be very wary of clicking the "Apply" button. "Apply" will apply the mapping you have chosen in the dialogue box to the pixels (assuming an 8-bit image or stack, or RGB stack).

(b) Sometimes it is useful to alter B&C automatically based on a section (a region of interest – ROI) rather than globally, note that the changes are applied across the whole image but it is the ROI that is used as the basis of this.

(c) For scientific images, changes to B&C must be made equally across the whole image.

### 2.2. Adding colour to a single grayscale image.

If you wish to add colour to a grayscale image so that you can open it as a colour image in other applications (such as Photoshop, Illustrator, PowerPoint etc) then you need to choose the colour and then save the image in RGB format.

2.2.1. Once the image is open select "Image" then "Lookup Table" then choose your colour. Note that the LUTs "Fire", "Ice" and "Spectrum" are *hot* meaning that colours are assigned based on pixel intensities – this might sometimes be useful. Once you have chosen the LUT the image will change accordingly. However, if saved in this state the colour information will be lost. In order to apply the colour information to the image the image type must be converted to RGB. Select "Image" then "Type" then "RGB Color". Next, save the image in a suitable format ("File", "Save As" and select type).

### N.B.

As a general rule images should always be saved as TIFFs. Never save as JPEG format since pixel information will be lost.

2.3. Working with channels – creating a merged image.

Often the image or stack contains more than one channel (colour) and after optimising B&C (Brightness and Contrast) of each separate channel you will want to create a merged or composite image.

- 2.3.1. If you opened the file as separate channels then you can create a merged image by choosing "Image", "Color", "Merge Channels".
- 2.3.2. If you opened the file without splitting the channels then you can first split the channels "Image", "Color" "Split Channels" or create a composite by "Image", "Color", "Make Composite".
- 2.3.3. Different LUTs can be applied to the channels before creating a composite, while if using "Merge Channels" then a fixed set of LUTs are selectable for each channel to be merged.
- 2.3.4. "B&C" can be modified for each split channel before merging, or for each channel within a image file than was not split by selecting the channel using the channel slider before altering "B&C".

The channel slider is marked with a "C"

2.4. Cropping an image.

Whilst it is desirable to zoom-in on the important information at the time of image capture, it is often necessary to crop images to highlight the important information, remove unimportant background, or to make images have the same proportions for reproduction.

2.4.1. Choose the appropriate crop tool from the menu (one of the first four buttons):

- 2.4.2. Select the area to be cropped and then select "Image" & "Crop". Note that once again this cannot be undone so it is best to crop a copy of the image ("Image" & "Duplicate").
- 2.4.3. To reuse exactly the same ROI (i.e. selection) on a second open image, click on the second image and then "Edit", "Selection" and "Restore Selection". This is very useful when you wish to crop a series of stacks to the exact same size.

### 2.5. Addition of a scale bar.

All micrographs should have a scale bar. This is typically very easy to apply in Fiji since the image is loaded with the appropriate metadata (i.e. the size of a pixel in micrometres). In some cases, depending on source of the image and possible bugs in the plugin used to read the file, the metadata may not be read correctly and the scale bar applied may be incorrect. Always look at your images and ask yourself if the scale bar looks correct!

2.5.1. To insert a scale bar first check if a scale has been attributed to the image from the metadata when the image was opened. Select "Analyze" and then "Set Scale", this dialogue box shows the number of pixels per micron as provided by the metadata. If no metadata is included with the file then the positions within the dialogue box will be empty. Check that the values are as you would expect given the imaging set up. For example, micrograph captured using a camera with pixel size of 6.5  $\mu$ m x 6.5  $\mu$ m and a 100X objective each pixel will be 0.065  $\mu$ m and therefore 1  $\mu$ m = 15.385 pixels. Once the information in the dialogue box is correct click "OK". Selecting "Global" will propagate this scale to all open images.

○○○ Set	Scale							
Distance in pixels:	15.385							
Known distance:	1.00							
Pixel aspect ratio:	1.0							
Unit of length:	micron							
Click to Remove Scale								
🗆 Global	Global							
Scale: 15.385 pixels/micron								
Help	Cancel OK							

2.5.2. Secondly, choose and stamp the scale bar. Select "Analyze" and then "Tools" followed by "Scale Bar" to show the dialogue box.

000 Scale	Bar Plus
Width in microns:	5
Height in pixels:	10
Font size:	14
Color:	White \$
Background:	None ‡
Location:	Lower Right \$
□ Bold Text ( □ Serif Font ( ☑ Label all slic	☑ Hide Text ☑ Overlay :es
(	Cancel OK

Here one chooses the thickness of the bar, its position, colour, presence/absence of text, and font style. Most journals prefer sans-serif fonts so this is default. Specific requirements will determine whether text is required. "Overlay" is recommended otherwise the bar will have the colour of the image. Equally, label all slices is recommended so that the scale is present throughout the image stack if appropriate. The location can be specified automatically in one of four positions, or using the position of the cursor. If you want to set by cursor position then the position should be marked using the point tool before opening the "Scale Bar" dialogue box:

### 2.6. Addition of a time stamp.

Time-lapse image stacks are best presented with a running clock so that viewers can see the real-time acquisition. There are two options to perform this: "Image", "Stacks", "Time Stamper" is one, but is not recommended. There is no preview function so getting the position and size correct is difficult and it is generally less flexible than the second option. The second option is "Image", "Stacks", "Series Labeler". Note that neither of these annotations can be undone (as often the case in Fiji). Therefore, all annotations should be performed on a copy (see paragraph 2 above – "Image" & "Duplicate"). Once the time-lapse stack is labelled it can be exported as a movie or individual images can be saved for presentation as a 2D figure composite.

2.6.1. Open the "Series labeller" dialogue: "Image", "Stacks", "Series Labeler". Choose the appropriate options within the dialogue box and enter the real time values in the text boxes. Select "Preview" at the bottom left of the dialogue box to assist with the size and placement of the time stamp. Generally the start time is "0", enter the interval between each image and generally choose to label every single image i.e. "Every n-th" is set to "1". "First" and "Last" are generally set automatically according to the size of the stack. Remember most journals prefer sans-serif fonts.

NB: If adding both a time stamp and a scale bar, add the time stamp first, otherwise the scale bar will be deleted, or flatten stack after adding scale bar ("Image", "Overlay", "Flatten").

000	Series Labeler
General Setting	5
Stack Type	time series or movie 🔹 🛊
Units Formattin	g
Label Format	Custom Format 💠 Decimal Places 0
Custom Suffix	sec
Start/Stop/Inter	val of Stack
Startup	0 First 0
Interval	4 Last 44
Every n-th	1
Location & Font	
X 20	Y 20 Location Presets Custom +
Font Settin	gs Font Color 🛛 Background
4: 173.85x173.85 microns (1024: 24 sec	1024) RGB 176M8 Werger-1
1. E.	
	* <i>1</i>

Example of a time series stack (of a photoconversion experiment using mitomEos) with time stamp in the top left and scale bar in the bottom right.

2.7. Addition of a z-position stamp.

Z-stacks to be presented as a movie are best labelled as described in paragraph 2.4. The only difference is that "z-stack" is chosen in the first drop-down menu

within the "Series Labeler" dialogue box. Again the labelled stack can be presented as a movie or as a series of 2D images in a composite figure.

2.8. Adding annotations to the image – using ROI manager.

Adding annotations (in this example, arrows) is achieved using the arrow tool and working with the ROI manager. Annotations are added to an overlay and are vector graphics as opposed to raster graphics and therefore are not affected by scaling, meaning they do not become pixelated. Only one ROI can exist at one time without adding them to the ROI manager.

2.8.1. To start it is best to make sure the arrow tool is added to the tool bar by right-clicking on the "line selection tool" which will add the arrow tool to the tool bar if it is not already in place. In the screen capture below, the arrow tool is in place.

000				(Fiji	ls Ju	st) In	nageJ							
IO I V	/ 4	+.	*	А	Q	ধ্য		Dev	Stk	LUŢ	Ø	b	\$ A	≫
(Fiji Is Just) Image 2.09-rc-2/1.49b; Java 1.6.0_65 [64-bit];														

2.8.2. Next, double click on the arrow tool to open the dialogue box:

ン <u>~</u> <u>∡</u> agel 2.0.0-r	(Fiji Is Just) Imag -+- <sup>×</sup> A	еј ] <u>Dev</u> Stk Lu [64-bit];	0 1 8 1
 000	Arrow Tool		
Width:	-	7	
Size:		10	
Color:	White ‡		
Style:	Filled \$		
□ 0 □ D ☑ K	utline ouble head eep after adding to o	overlay	
	Cancel	ОК	

2.8.3. Now open the ROI manager – "Analyse", "Tools", "ROI Manager". Be sure to select "Show All" in the dialogue box otherwise you will not be able to see each arrow you draw: if "Show All" is not selected you will only see the last arrow added. I also find it useful to select "Labels" to add a number to each arrow while they are being edited.

000 ROI Man	ager
	Add [t]
	Update
	Delete
	Rename
	Measure
	Deselect
	Properties
	Flatten [F]
	More »
	Show All
	<b>✓</b> Labels

2.8.4. Draw the arrow you want, and arrange it as required. Now click on the "Add [t]" button, or simply press "t" on the keyboard.

2.8.5. If you want to add more than one arrow to the same slice (image) then you will have to add the ROI to the overlay by pressing "b" on the keyboard.

2.8.6. If you are working with a stack (time series or z) and want to label more slices with the same arrow then you must work with each arrow in turn on the different slices i.e. do not press "b" to add the arrow to the overlay (see bewlow) otherwise it will no longer be an ROI and will not be transferable

to another slice. Instead, draw the arrow, position it (it can be positioned as long as it is indicated by little dots: when the arrow is positioned correctly press "t" and then move to the next slice using the ">" key (N.B. This is <u>not</u> an arrow key). The arrow can be re-positioned if necessary and added to the "ROI Manager" by pressing "t". All ROIs are still editable and can be deleted if necessary.

- 2.8.7. When you are happy with all the annotations then you must add the ROIs as an overlay this is done by selecting each ROI in "ROI Manager" and typing "b". At this stage, if you save the stack then the overlay will be saved with it. <u>However</u>, the overlay is only seen by ImageJ, not by other programs, i.e. if you open in Photoshop then the arrows will not be present.
- 2.8.8. To imbed the arrows (overlay) in the image (and here it is final, so work on a copy "Image" & "Duplicate"), then you must click on "Flatten" in the "ROI Manager" or type "f".

### 2.9. Exporting as a movie file.

Stacks are easily saved as movies in Fiji: "File", "Save As", "AVI". Uncompressed is recommended but this depends on the downstream use. As a Mac user I always convert my movies to QuickTime movies so any compression is performed at that stage. The only other option is to choose a playback rate (in frames per second). The beauty of having already added the time stamp is that this playback rate is less important (i.e. does not need to be provided to understand the movie) since viewers can see the real-time acquisition rate stamped on the movie. Choosing playback rate may require some trial and error.

### 3. Analysis of mitochondrial dynamics – qualitative and quantitative.

The capture of good quality images to show mitochondrial morphology and dynamics is itself a challenge but beyond this the community is becoming more demanding in requesting quantification of morphology and dynamics. Some quantification can be relatively simple to achieve i.e. numbers of mitochondria, numbers of long versus round mitochondria, number of mitochondrial clusters versus number of single mitochondria. Some of these can be obtained by simple counting, but it is best if reproducible non-biased criteria can be applied. ImageJ can be used for these measurements but there is no time to go into all possibilities. Here two quantitative assays are presented – co-localisation to quantify the extent of mitochondrial movement; and particle counting to measure the numbers, sizes and morphology indicators of mitochondria.

### 3.1. Co-localisation for qualitative and quantitative representation of dynamics.

3.1.1. Qualitative representation: in it's simplest qualitative form, one can capture two images a precise time frame apart, apply different LUTs to each image and then merge them. This method gives a visually very striking representation of mitochondrial motility.

In the figure above two images were captured, 1 second apart, of mitochondria moving in an Arabidopsis leaf cell. The first image (A) was false-coloured green and the second (B) false-coloured red, using LUTs (see section 2.2). The resulting merged image is shown in (C) where mitochondria that moved within the 1 s period are coloured green or red while those in the same position are yellow.

- 3.1.2. In order to compare treatments or conditions that may affect motility it is beneficial to be able to generate quantitative data. One relatively simple method to do this is to use co-localisation algorithms within ImageJ, for example "Analyze", "Colocalization", "Colocalization Threshold". This method provides colocalisation statistics and an image similar to that above but with colocalised pixels in white. Although one of the biggest problems with co-localisation is thresholding of the images this is not relevant here since we are using two images captured with exactly the same parameters and not the more usual colocalisation of one fluorophore with another that would typically be captured with different settings (illumination and detector).
- 3.1.3. Open the two images to be compared, i.e. images captured at two defined time points, or extract these single time points from a stack (this can be achieved easily using "Image", "Stacks", "Tools", "Slice Keeper"

and defining the slice to be copied from the stack as a single image – do this for each of the two images to be used for colocalisation analysis).

3.1.4. Perform colocalisation "Analyze", "Colocalization", "Colocalization Threshold" - the following dialogue box will appear:

○ ○ ○ Colocalisation Thresholds			
Channel 1	5.tif \$		
Channel 2	6.tif \$		
Use ROI	None ‡		
Channel Combination	Red : Green ‡		
<ul> <li>✓ Show Colocalized</li> <li>✓ Use constant integer</li> <li>✓ Show Scatter plot</li> <li>✓ Include zero-zer</li> <li>✓ Set options</li> <li>See: http://uhnrese</li> </ul>	d Pixel Map ensity for colocalized pixels t o pixels in threshold calculation arch.ca/wcif/imagej Cancel OK		

3.1.5. Select the two images as "Channel 1" and "Channel 2" and select all other options, for example, "Set options" must be enabled to view the statistics. Selecting "Use constant intensity for colocalised pixels" means that they will appear white in the "Colocalised Pixel Map".

Above, the two images to be subject to colocalisation analysis, note that the images do not need to be coloured, grayscale works fine. Below, the "Colocalised Pixel Map":

000					Results				
Images	Mask	ZeroZero	Rtotal	R <threshold< th=""><th>M1</th><th>M2</th><th>%Ch1 Int &gt; thresh</th><th>%Ch2 Int &gt;thresh</th><th></th></threshold<>	M1	M2	%Ch1 Int > thresh	%Ch2 Int >thresh	
5.tif & 6.tif	ROI0	incl.	0.519	0.006	0.7674	0.7742	72.67%	73.39%	ľ

3.1.6. A variety of statistics can be generated and output to a results sheet. See <a href="http://www.uhnresearch.ca/facilities/wcif/imagej/colour\_analysis.htm">http://www.uhnresearch.ca/facilities/wcif/imagej/colour\_analysis.htm</a> for a description of these. For the purposes of quantifying mitochondrial motility we can quote Pearson's or Mander's coefficient and the "% Intensity above threshold colocalized". This later value, like the Mander's coefficient, is presented for each channel and is likely to be different in the two cases since the two images (channels = time points in this case) are of course different. To provide a value for colocalisation of mitochondria as an indicator of lack of movement an average of the two channel values can be used. For the example image, showing little movement, this value is 73.03%.

- 3.2. Counting and measuring particles (mitochondria).
  - 3.2.1. To simply count particles in an image one can use the "Process" "Find Maxima" command. The settings require optimisation of the "Noise Tolerance" value for each image simply vary this value while previewing point selection (i.e. have "Preview point selection" checked).

O O O Find Maxima	
Noise tolerance: 100	
Output type: List +	
<ul> <li>✓ Exclude edge maxima</li> <li>□ Light background</li> <li>✓ Preview point selection</li> </ul>	
175 Maxima	
Help Cancel OK	
o o putative_mutant_time_series_2.czi kept st. 1040x856 pixels; 16-bit; 1.7MB	ick-1
¢ \$	
6 6 0° 0°	•• <sup>•</sup> •• •
6 •	·
* 200	· · · · * *
о <sup>р</sup>	
• • • • • • • • • • • • • • • • • • •	<u>, 0</u>
° ° ° °	· · · ·
• • •	*
+ ? * +	+ ```a +
	0 ÷
	•
	÷. •

- 3.2.2. Selecting "Count" from the "Output type" dropdown menu returns the number of maxima found. Selecting "List" returns the x and y coordinates of each maximum. "Find Maxima" does not work on stacks but the "Find Stack Maxima" macro can be used.
- 3.2.3. Select "Plugins", "Macros", "Record" and paste the macro below (also available here: <u>http://imagej.nih.gov/ij/macros/FindStackMaxima.txt</u> into the window then click on "Create" to open the editor. Then "File" "Save As" to save the macro. Then select "Plugins" "Macros" "Install" to find and install the macro file you created. This macro should then appear when you select "Plugins", "Macros".

```
// Find Stack Maxima
11
// This macro runs the Process>Binary>Find Maxima
// command on all the images in a stack.
 Dialog.create("Find Maxima");
 Dialog.addNumber("Noise Tolerance:", 5);
 Dialog.addChoice("Output Type:", newArray("Single
Points", "Maxima Within Tolerance", "Segmented Particles",
"Count"));
  Dialog.addCheckbox("Exclude Edge Maxima", false);
 Dialog.addCheckbox("Light Background", false);
 Dialog.show();
 tolerance = Dialog.getNumber();
 type = Dialog.getChoice();
 exclude = Dialog.getCheckbox();
  light = Dialog.getCheckbox();
 options = "";
  if (exclude) options = options + " exclude";
  if (light) options = options + " light";
  setBatchMode(true);
  input = getImageID();
 n = nSlices();
  for (i=1; i<=n; i++) {</pre>
     showProgress(i, n);
     selectImage(input);
     setSlice(i);
     run("Find Maxima...", "noise="+ tolerance +"
output=["+type+"]"+options);
     if (i==1)
        output = getImageID();
    else if (type!="Count") {
       run("Select All");
       run("Copy");
       close();
       selectImage(output);
       run("Add Slice");
       run("Paste");
    }
  }
  run("Select None");
  setBatchMode(false);
```

- 3.2.4. To count and measure particles one can use "Analyze", "Analyze Particles".
- 3.2.5. First, the image must be subjected to threshold using "Image" "Adjust" "Threshold". Thresholding is not easy to perform and almost every image

is different, at least when imaging plant tissue. I imagine thinner, flatter animal cells in culture are much easier to threshold. As ever, the ease of thresholding and success of the method depends on the cleanliness and contrast of the starting image. Images with uniformly bright objects against a uniformly dark background are the easiest to threshold. Such images are not easy to capture! It is unlikely all objects (mitochondria) will be included after thresholding but the analysis of multiple images and averaging of the results should overcome this problem. Equally, thresholding smaller ROIs within the image may work better if trying to threshold large bright objects and small dim objects in the same image.

3.2.6. Threshold image: "Image", "Adjust", "Threshold". In this example red is chosen to highlight the particles to be retained. Once the appropriate threshold is set, click "Apply". It is rare that "Auto" results in an acceptable threshold. Clicking "Set" allows you to input values for the upper and lower thresholds rather than using the sliders.

000	Threshold			
0.90 %		631 1263		
Intermodes	‡ Red	\$		
🗹 Dark background 🛛 Stack histogram				
Auto	Apply Reset	Set		

It is often helpful to "Process", "Smooth" the image before thresholding to reduce the chances that small speckled areas will be identified:

Original

Smoothed

3.2.7. N	ow you can	"Analyze",	"Analyze	Particles":
----------	------------	------------	----------	-------------

○ ○ ○ Analyze Particles				
Size (µm^2): 0-Infinity Pixel units				
Circularity: 0.00-1.00				
Show: Outlines ‡				
<ul> <li>Display results</li> <li>Exclude on edges</li> <li>Clear results</li> <li>Include holes</li> <li>Summarize</li> <li>Record starts</li> <li>Add to Manager</li> <li>In situ Show</li> </ul>				
Help Cancel OK				

If the image has been calibrated for size then the size of the particles to be included in the analysis can be defined in the first text box. This is useful to exclude wrongly identified due to thresholding or to exclude clumps of particles that cannot be resolved as individual particles.

"Circularity" refers to how round the particles are – this value will be reported in the "Results" window ("Window", "Results") and is useful for assigning a quantitative value to the mitochondrial morphology. Particles are identified with an ellipse; a value of 1 is an ellipse that is a perfect circle, while a value of 0 is a straight line. If you were dealing with uniform particles this could be a means to filter, for mitochondria it is instead a useful measure of morphology.

The results to be reported are selected using the "Analyze", "Set Measurements" dialogue box:

000	Set Measurements
<ul> <li>Area</li> <li>Standard deviation</li> <li>Min &amp; max gray value</li> <li>Center of mass</li> <li>Bounding rectangle</li> <li>Shape descriptors</li> <li>Integrated density</li> <li>Skewness</li> <li>Area fraction</li> </ul>	<ul> <li>Mean gray value</li> <li>Modal gray value</li> <li>Centroid</li> <li>Perimeter</li> <li>Fit ellipse</li> <li>Feret's diameter</li> <li>Median</li> <li>Kurtosis</li> <li>Stack position</li> </ul>
Limit to threshold Invert Y coordinates Add to overlay	Display label Scientific notation
Redirect to:	lone 🗘
Decimal places (0-9): 3	
	Help Cancel OK

The following is taken from the ImageJ user guide

http://rsbweb.nih.gov/ij/docs/guide/146.html

*«* **Area** Area of selection in square pixels or in calibrated square units (e.g., mm2,  $\mu$ m2, etc.) if Analyze  $\triangleright$  Set Scale... $\downarrow$  was used to spatially calibrate the image.

**Mean gray value** Average gray value within the selection. This is the sum of the gray values of all the pixels in the selection divided by the number of pixels. Reported in calibrated units (e.g., optical density) if Analyze  $\triangleright$  Calibrate...↓ was used to calibrate the image. For RGB images, the mean is calculated by converting each pixel to grayscale using the formula gray = (red + green + blue)/3 or gray = 0.299 × red + 0.587 × green + 0.114 × blue if Weighted RGB Conversions is checked in Edit  $\triangleright$  Options  $\triangleright$  Conversions...↑

*Standard deviation* Standard deviation of the gray values used to generate the mean gray value. Uses the Results Table↑ heading **StdDev**.

*Modal gray value* Most frequently occurring gray value within the selection. Corresponds to the highest peak in the histogram. Uses the heading **Mode**.

Min & max gray level Minimum and maximum gray values within the selection.

*Centroid* The center point of the selection. This is the average of the x and y coordinates of all of the pixels in the image or selection. Uses the **X** and **Y** headings.

*Center of mass* This is the brightness-weighted average of the x and y coordinates all pixels in the image or selection. Uses the **XM** and **YM** headings. These coordinates are the first order spatial moments.

**Perimeter** The length of the outside boundary of the selection. Uses the heading **Perim.** With IJ 1.44f and later, the perimeter of a composite selection is calculated by decomposing it into individual selections. Note that the composite perimeter and the sum of the individual perimeters may be different due to use of different calculation methods.

**Bounding rectangle** The smallest rectangle enclosing the selection. Uses the headings **BX**, **BY**, **Width** and **Height**, where **BX** and **BY** are the coordinates of the upper left corner of the rectangle.

**Fit ellipse** Fits an ellipse to the selection. Uses the headings **Major**, **Minor** and **Angle**. **Major** and **Minor** are the primary and secondary axis of the best fitting ellipse. **Angle** is the angle between the primary axis and a line parallel to the X-axis of the image. The coordinates of the center of the ellipse are displayed as **X** and **Y** if Centroid is checked. Note that ImageJ cannot calculate the major and minor axis lengths if *Pixel Aspect Ratio* in the Analyze  $\triangleright$  Set Scale...↓ dialog is not 1.0. There are several ways to view the fitted ellipse:

- 1. The Edit ▷ Selection ▷ Fit Ellipse↑ command replaces an area selection with the best fit ellipse.
- 2. The DrawEllipse macro draws (destructively) the best fit ellipse and the major and minor axis.
- Select *Ellipses* from the *Show:* drop-down menu in the particle analyzer (Analyze ▷ Analyze Particles...↑) and it will draw the ellipse for each particle in a separate window.

Shape descriptors Calculates and displays the following shape descriptors:

Circularity 4π × [Area][Perimeter]2with a value of 1.0 indicating a perfect circle. As the value approaches 0.0, it indicates an increasingly elongated shape. Values may not be valid for very small particles. Uses the heading Circ. Aspect ratio The aspect ratio of the particle's fitted ellipse, i.e., [Major Axis][Minor Axis]. If Fit Ellipse is selected the Major and Minor axis are displayed. Uses the heading AR. Roundness 4 × [Area]π× [Major axis]2or the inverse of Aspect Ratio. Uses the heading Round. Solidity [Area][Convex area]; Note that the Edit ▷ Selection ▷ Convex Hull↑ command makes an area selection convex.

Feret's diameter The longest distance between any two points along the selection boundary, also known as maximum caliper. Uses the heading
Feret. The angle (0--180 degrees) of the Feret's diameter is displayed as
FeretAngle, as well as the minimum caliper diameter (MinFeret). The starting coordinates of the Feret diameter (FeretX and FeretY) are also displayed (*see also* Feret's Diameter macro and Chamfer distances and Geodesic diameters plugin).
Integrated density The sum of the values of the pixels in the image or selection. This is equivalent to the product of Area and Mean Gray Value. With IJ 1.44c and later, Raw integrated density (sum of pixel values) is displayed under the heading RawIntDen when Integrated density is enabled. The Dot Blot Analysis tutorial demonstrates how to use this option to analyze a dot blot assay.

*Median* The median value of the pixels in the image or selection.

*Skewness* The third order moment about the mean. The documentation for the Moment Calculator plugin explains how to interpret spatial moments. Uses the heading **Skew**.

*Kurtosis* The fourth order moment about the mean. Uses the heading Kurt.

**Area fraction** For thresholded images is the percentage of pixels in the image or selection that have been highlighted in red using  $Image \triangleright Adjust \triangleright Threshold... [T]\uparrow$ . For non-thresholded images is the percentage of non-zero pixels. Uses the heading **%Area**.

*Stack position* The position (slice, channel and frame) in the stack or hyperstack of the selection. Uses the headings **Slice**, **Ch** and **Frame**.

N.B.: For line selections the heading **Length** is created. For straight line selections, **Angle** is recorded even if *Fit Ellipse* is unchecked. Also, note that measurements that do not apply to certain selection types may be listed as *NaN*, *Infinity* or -*Infinity*. »

3.2.8. By selecting "Add to Manager" the identified particles are added to "Analyze", "Tools", "ROI Manager". Selecting one of the labelled particles in "ROI Manager" will highlight the corresponding particle in the overlay image. Mis-identified particles can be deleted using this dialogue box.

3.2.9. To make a composite image of the original and the identified particles outline image, copy the outline, "Edit", "Copy", then "Edit", "Paste Control" and select "Add" from the pulldown menu:

000 P	aste Control	
Transfer Mode:	Add	\$

Then click on the original image and "Edit", "Paste". Finally, select "Image", "Overlay", "Show Overlay".

Original

Threshold